

Health o meter[®]
Professional **MECHANICAL CHAIR SCALE**

- Steel base for strength and durability
- Mechanically measures weight
- Folding foot rest and side arms for easy patient access and comfort
- Heavy-duty wheels for mobility
- 2 Year Warranty
- 440lb/220kg weight capacity

654-445KLEach



Health o meter[®]
Professional **DIGITAL CHAIR SCALE**

- Steel base for strength and durability
- Display: Large 1.25" backlit LCD display
- LB / KG Conversion, Zero Out / Tare, Auto Zero, Auto Off
- Powered by a 6V rechargeable battery and battery charger included
- Folding foot rest and side arms for easy patient access and comfort
- Heavy-duty wheels for mobility
- 2 Year Warranty
- 600lb/270kg weight capacity

654-594KLEach



Health o meter[®]
Professional **DIGITAL STAND ON SCALE WITH HANDRAILS**

- LB / KG Conversion, LB /KG Lock Out, BMI, Tare, Hold/Release, Reweigh, Auto Zero, Auto Off
- Display: LCD, 180° Swivel
- Oversized 26"X 22" Platform
- Patient can hold onto handrails for stability while scale accurately calculates weight
- Wheels for ease of movement
- Powered by 6 AA-cell batteries (not included) or 120V AC adapter (included)
- 2 Year Warranty
- 1000lb/454kg weight capacity
- Optional height rod available. Part # 654-STROD

654-2101KL .. ScaleEach



Health o meter[®]
Professional **DIGITAL WHEELCHAIR RAMP SCALE**

- LB / KG Conversion, LB /KG Lock Out, BMI, Tare, Hold/Release, Reweigh, Auto Zero, Auto Off
- Display: LCD, 180° Swivel
- Platform size 29.25" x 28.25"
- Versatile scale can be used with a regular chair, wheelchair, or as a stand on scale with support handrail
- 2 wheels for ease of movement. Ramp does not fold.
- Powered by 6 AA-cell batteries (not included) or 120V AC adapter (included)
- 2 Year Warranty
- 1000lb/454kg weight capacity



654-2500KLEach

Health o meter[®]
Professional **DIGITAL FOLDING WHEELCHAIR RAMP SCALE**

- LB / KG Conversion, LB /KG Lock Out, BMI, Tare, Hold/Release, Reweigh, Auto Zero, Auto Off
- Display: LCD, 180° Swivel
- Platform size 32"x36" accommodates standard and bariatric wheelchairs
- Versatile scale can be used with a regular chair, wheelchair, or as a stand on scale with support handrail
- Scale platform folds for convenient storage with 4 wheels for ease of movement
- Powered by 6 AA-cell batteries (not included) or 120V AC adapter (included)
- 2 Year Warranty
- 1000lb/454kg weight capacity



654-2600KL Single side rampEach
654-2610KL Dual side rampEach

Health o meter[®]
Professional **DIGITAL PORTABLE WHEELCHAIR SCALE**

- LB / KG Conversion, LB /KG Lock Out, BMI, Tare, Hold/Release, Reweigh, Auto Zero, Auto Off
- Scale folds and includes built-in handle for easy portability
- Platform expands from 25-1/2" to 50" to accommodate bariatric and standard wheelchairs
- Powered by 6 AA-cell batteries (not included) or 120V AC adapter (included)
- 800lb weight capacity
- 2 Year Limited Warranty



654-2400KLEach

Health o meter[®]
Professional LARGE PLATFORM SCALE 8" DIAL

Mechanical personal scale with large dial and non-slip platform.
 Capacity 400lbs.

654-160KL Each



Health o meter[®]
Professional BALANCE BEAM WEIGH SCALE

- Mechanically measures weight in pounds and kilograms.
- Mechanical rotating poise bars allow for easy reading of weight in lb or kg
- Engraved numerals on the poise bars can withstand frequent use, providing long-lasting readability
- Easy-to-use height rod provides height measurements from 23" to 84"/60cm-213cm
- 500 lb / 200 kg weight capacity - with included counterweights
- 10 Year Limited Warranty

654-450KL Each



Health o meter[®]
Professional DIGITAL EYE LEVEL SCALE

The Elevate EMRscale digitally measures height, weight and BMI and transmits those measurements via its Connectivity Kit. This Connectivity Kit is an external adapter that allows the scale to interface with Welch Allyn Connex Vital Signs Monitors, Connex[®] Spot Monitors, and Connex[®] Integrated Wall Systems. From the Welch Allyn Connex[®] device, the data can then be displayed, saved, and transferred into EMR systems. Additional Features and Benefits of the Elevate EMRscale™

- Fast and simple digital calculations of weight, height and BMI speeds process of vital sign measurement
- Digital transmission of vital signs reduces the risks of error that can occur with manual data entry
- Low Profile Platform requires only minimal step up, offering patient a more comfortable and stable weighing experience
- High 600lb capacity

654-ELEVATE-C Each



DIGITAL CHAIR SCALE

- This scale is designed for the mobility-challenged who have a difficult time standing on a conventional scale.
- Equipped with a plastic seat and large 1" LCD display.
- Swing-up arms enable easy access and exit.
- Movement compensation technology eliminates inaccuracies caused by involuntary movement during the weighing process.
- Units of measure in lb only, kg only and lb/kg.
- Includes internal rechargeable battery, low battery indicator (auto shut off), BMI (body mass index) function, reweigh, hold function and tare function. AC Adapter is optional #MDR133077

320-MDR119114Each
 320-MDR133077.....Each



DIGITAL WHEELCHAIR SCALE WITH RAMP

- Non-skid platform
- Movement compensation technology
- LCD display
- Functions: BMI, Hold, and LB, KG and LB/KG measurement
- EMR capable, USB
- Powered by 6AA batteries (included)
- Optional AC adapter available
- Weight capacity: 1,000 lbs. (454 kg)

320-MDR141447.....Each



BALANCE BEAM MECHANICAL WEIGH SCALE

- Mechanical physician scale features rear wheels for easy portability
- Durable, accurate and easy to operate, move and store
- Weight can be read from either the front or back
- Slip-resistant base cover is removable for easy cleaning
- Height rod: 24" - 84" (60 - 212 cm)
- Graduation: 0.25 lb (0.1 kg)
- Unit of measure: lbs/kg
- 2-year warranty

320-MDR102613.....Each



seca DIGITAL FLOOR SCALE

The classic, clean line of the basic model turns heads and the high-quality rubber coating feels good to the touch. It's also pleasantly easy to read the weights in the display because of the large LCD figures.

- Capacity: 330 lbs / 150 kg
- Graduation Weight: 0.2 lbs / 100 g
- Weight: 3.3 lbs / 1.5 kg

407-803 Each



seca DIGITAL COLUMN SCALE

With its durable cast iron base and high weight capacity of 660 pounds, the seca 703 is suitable for all patient types. The low level platform is slip resistant giving patients a firm foothold for extra stability and safety. The integrated measuring rod seca 220 turns the scale into a complete weighing and measuring station with automatic BMI calculation. The practical transport castors permit mobile use for convenient weighing in various locations. With everything primed for the future and integration in any EMR system, measurements can be transmitted wirelessly via the seca 360° wireless system.

- Heavy duty cast iron base
- Slip resistant platform with an anti-tip column design
- BMI, Pre-TARE, HOLD, Auto-HOLD, lbs/kg switch-over, damping, SEND/PRINT
- Included measuring rod seca 220
- At the press of a button, measurements can be wirelessly transmitted to a seca 360° Wireless digital printer, PC or EMR system

407-703 Each

