

↑ EZee Life SHOWER & COMMODE CHAIR - STANDARD



798-180 Each

- Adjustable height
- Light weight aluminum NO RUST
- Slide on pail, removable
- Footrest, adjustable height & removable
- Flip back arms
- Anti tipper bars on rear
- Calf strap included
- Four 5" locking casters
- Super soft seat (removable insert)
- Seat belt
- Optional 24" rear wheels are available
- All adjustments are tool-less
- 2 year warranty

Specifications	
Weight Capacity	300 lbs (136.4 kgs)
Rear Wheels	5" standard or optional 24"
Front Casters	Locking casters
Seat Width	18"
Seat Depth	17"
Seat to Floor Height	19¼" to 22¼"
Back Height	18"
Arm Length	11" (from seat)
Product Weight	30 lbs (13.6 kgs)
Rear Internal Width	15¾"
External Width	19"
External Length	37" to 38 ½"
Clearance (Floor to Rail)	17" to 20"

↑ EZee Life SHOWER & COMMODE CHAIR - BARIATRIC



798-185 .. 20" .. Each

798-188 .. 22" .. Each

- Adjustable height
- Light weight aluminum NO RUST
- Slide on pail, removable
- Footrest, adjustable height & removable
- Flip back arms
- Anti tipper bars on rear
- Calf strap included
- Four 5" locking casters
- Super soft seat (removable insert)
- Seat belt
- Optional 24" rear wheels are available
- All adjustments are tool-less
- Available in 20" width or 22"Width
- 2 year warranty

Specifications	
Weight Capacity	400 lbs (189 kgs)
Rear Wheels	5" standard or optional 24"
Front Casters	Locking casters
Seat Width	20" or 22"
Seat Depth	17"
Seat to Floor Height	19¼" to 22¼"
Back Height	18"
Arm Length	11" (from seat)
Product Weight	30 lbs (13.6 kgs)
Rear Internal Width	17¾" Or 19¾"
External Width	21" or 23"
External Length	37" to 38 ½"
Clearance (Floor to Rail)	17" to 20"

↑ EZee Life SHOWER & COMMODE CHAIR - TILT



798-190 .. 18" .. Each

798-195 .. 20" .. Each

- Adjustable height
- Light weight aluminum NO RUST
- Slide on pail, removable
- Footrest, adjustable height & removable
- Flip back arms
- Anti tipper bars on rear
- Calf strap included
- Four 5" locking casters
- Super soft seat (removable insert)
- Seat belt
- Optional 24" rear wheels are available
- Adj. headrest up/down/forward/back
- Tilt Range 0-20° Posterior & 0-10° Anterior
- All adjustments are tool-less
- Available in 18" & 20" widths
- 2 year warranty

Specifications	
Weight Capacity	300 lbs (136.4 kgs)
Rear Wheels	5" standard or optional 24"
Front Casters	Locking casters
Seat Width	18" or 20"
Seat Depth	17"
Seat to Floor Height	21¼" to 24¼"
Back Height	18"
Arm Length	11" (from seat)
Product Weight	36 lbs (16.4 kgs)
Tilt Range Posterior	0-20°
Tilt Range Anterior	0-10°
Rear Internal Width	15¾" or 17¾"
External Width	19" or 21"
External Length	37" to 38 ½"
Clearance (Floor to Rail)	19" to 22"

PVC SHOWER/COMMUNE CHAIR WITH FIXED ARMS

- 18" Internal width
- Thermoformed seat for easy cleaning (220 model)
- 7 Quart pail with lid
- 3" Tente total lock casters on 220 model, 3" standard lock casters on 200 model
- Available in Dark Blue (B), Dark Green (G), or Mauve (M) mesh backrest

798-200(colour)	Standard Shower/Commune Chair	Each
798-220(colour)	Deluxe Shower/Commune Chair	Each



PVC SHOWER/COMMUNE CHAIR WITH DROP-ARM

- 18" Internal Width
- Double Drop-arm
- Thermoformed seat for easy cleaning
- 7 Quart pail with lid
- 3" Tente total lock casters
- Available in Dark Blue (B), Dark Green (G), or Mauve (M) mesh backrest

798-250(colour)	Shower/Commune Chair	Each
---------------------------	--------------------------------	------



PVC BARIATRIC SHOWER/COMMUNE CHAIRS

- Thermoformed seat for easy cleaning
- 7 QT Pail
- Fixed Curved Arms
- 3" or 5" Tente Casters, 4 Locking
- Weight Capacity: 350 lbs or 400 lbs
- 21" or 24" Seat Widths

798-260(colour)	Fixed Curved Arms, 350 lb Weight Capacity, 21" Seat, 3" Casters
798-260-T5(colour)	Fixed Curved Arms, 350 lb Weight Capacity, 21" Seat, 5" Casters
798-270(colour)	Fixed Curved Arms, 400 lb Weight Capacity, 24" Seat, 3" Casters
798-270-T5(colour)	Fixed Curved Arms, 400 lb Weight Capacity, 24" Seat, 5" Casters



REPLACEMENT PARTS

798-C3TSB	3" Tente Caster with lock	Each
798-C3TS	3" Tente Caster no lock	Each
798-VACS-18	Thermoform Seat 18"	Each
798-VACS-21	Thermoform Seat 21"	Each
798-B	Seat Belt	Each
798-C7	7 Quart Commode Pot	Each



DELUXE DROP ARM WHEELED COMMUNE

Commune is a truly versatile commode that features 4 swivel casters, removable swing-away footrests, and a removable padded seat. What's more, patient transfers are fast and easy with the commode's fold-down padded armrests and front and rear locking casters. (6.6 L pail with cover and handle included)

A41-11120SV-1SF Each



AQUATEC OCEAN® COMMUNE/SHOWER CHAIR

The multi-functional Ocean can be used as a height adjustable commode chair, shower chair or bedside commode. Features include stainless steel construction, tension adjustable backrest, removable height adjustable swing-away footrests with adjustable heel loops, commode pan with lid, four 5" casters with wheel & swivel locking mechanism, ergonomically designed flip-up armrests. 300lb Capacity.

684-OCEAN 18" Width Each
 684-OCEANXL 21" Width Each



TILT SHOWER CHAIR/COMMUNE

- Height Adjustable Base
- Nickel Plated Epoxy Coated Frame
- 4 x 5" Rust Resistant Nylon Casters, 3 Total Lock
- Plastic Seat with Bucket and Lid
- Waterproof Mesh Back
- Flip Back Cantilever Arms with Water Proof Pads
- Gas Struts for Tilting
- Customize and Upgrade Options Available.

678-TILTCOMC 18" Seat Width Each
 678-TILTCOMC/20 20" Seat Width Each



ACTIVEAID® SHOWER CHAIR/ COMMUNE

- 250 lb weight capacity
- Seamless Ensolite™ padded seat and back
- Swingaway, removable arms
- Adjustable, removable footrests
- Push handle back
- Removable commode pan holder and pan
- Corrosion-resistant stainless steel construction
- Two dual-lock casters



649-6055S Each

ACTIVEAID® TILT-IN SPACE SHOWER CHAIR/COMMUNE

- Stainless steel frame with non-corrosive composite components
- Gravity Assist Tilting system (patented)
- Removable, height adjustable footrests
- Padded seat with anti-drip inderring
- Flip-back padded removable arm with locking system for weight shift
- Dual locking 5" Stainless steel/ composite casters
- 350lb capacity



649-2855S 18" Seat Width Each

649-2855S20 20" Seat Width Each

ACTIVEAID® BARIATRIC SHOWER CHAIR/ COMMUNE

- 600 lb capacity
- Stainless steel, high strength rigid frame
- Seamless Ensolite™ padded seat and back cushion
- Swingaway, removable padded arms
- Flip-up, removable foot platform
- Removable commode pan
- 20" rear wheels
- Toggle wheel locks
- 6" front casters
- Push handle back



649-720 20" seat width Each

649-722 22" seat width Each

649-724 24" seat width Each

649-726 26" seat width Each

ECONOMY SHOWER CHAIR/COMMODE

This shower chair can also double as a commode seat. It features a satin finished aluminum frame, full length plastic armrests, comfortable oversized padded toilet seat, rust-resistant 3" swivel casters, diagonal locking casters, and a comfortable backrest with cover (included).

A41-11114KD-1 Each



ECONOMY WHEELED COMMODE

Four swivel casters with diagonal locking brakes, angled back support, and full length plastic armrests makes this commode ideal for institutional or long term home users. This commode also features satin finished aluminum frame, plastic seat and cover, 3" rubber swivel casters, and a 7 litre plastic pail and cover. The back support can be knocked down to minimize storage space. As well, with the pail removed, it can be used as a toilet safety frame.

A41-11112-2 Each



STATIONARY COMMODE, ADJUSTABLE HEIGHT

- Deep snap on seat with lid.
- Powder coated steel welded construction increases strength and durability.
- Removable tool free back
- Plastic arms for added comfort
- Comes complete with a 12qt. commode bucket with carry handle, cover and splash shield.

320-MDS89664 Standard Width Each

320-G30216B X-Wide Width Each



Standard

X-Wide

BEDPAN & COMMODE PAIL ABSORBENT BAG

The bedpan cover is to be used in conjunction with a standard bedpan or over a standard commode pail. The cover has an absorbent pad that absorbs up to 500cc in organic liquid. The system provides an effective way to manage the collection, transportation, smell and disposal of organic waste. Bags are biodegradable and can be used over commode pails for the same purpose. 40cm x 52cm

320-MDS89664LINER Hygienic Bedpan & Commode Pail Cover 72/Box

