

## Evacusafe™ EXCEL EVACUATION CHAIR

Evacusafe Excel Evacuation Chairs are designed to enable the safe and comfortable evacuation of mobility impaired people in the event of an emergency, such as fire. The chairs are designed for ease-of-use on internal stairways and fire escapes as well as along flat surfaces.

Evacuation Chairs are suitable for a single person operation, even with up to 400 pounds on board. Mainly operated on a straight continuous run, they can also be used (with some training) on a mild spiral, though this will be best suited for a 2 person operation.

- Lightweight yet incredibly strong
- Can transport persons of up to 180Kg
- Handle can fold away and the chair can then be used as a temporary wheelchair
- Supportive FLAT seat, foot and armrests and padded headrest, the chair is very comfortable to stay in for extended periods
- Folds flat and can be stored on the wall with the hooks supplied along with an 'Evacuation Chair' wall sign
- No Lifting or manual handling required
- Tough, long lasting cover helps keep the equipment clean, protected and ready for use
- Easy to set up and use
- Online training is available via the manufacturer's online training site for a nominal fee
- Ideal for downward stairway evacuation
- Rubberized friction activated tracks
- 4 wheels for optimum manoeuvrability
- Foot operated brake
- Built in carry handles
- Quick release safety harness
- Removable covers

### Specifications:

Weight ..... 13.6kg  
 Load Capacity ..... 182kg  
 Height (folded/open) . . . 1130mm / 1320mm  
 Depth (folded/open) . . . 290mm / 990mm  
 Width ..... 650mm  
 Length in use. .... 1430mm

C57-EXCEL ..... Each

Swift,  
Smooth,  
Safe



### THREE POSITION RECLINER WITH TRAY

- 3 positions: upright, deep recline and elevated legrest/footrest.
- Durable, fully enclosed side panels securely attached to protect clothing and skin.
- Adjustable seat straps provide maximum comfort.
- Self adjusting suspension in the seat prevents bottoming out.
- Large, blow moulded tray can be stored on the side of the chair.
- 5" casters (2 with locks) make it easy to move and secure the recliner.
- Wide, comfortable armrests are angled to provide support in any reclined position.
- Neutral coloured, sturdy steel frame is easy to clean and maintain.
- Moisture barrier on seat prevents seepage.
- Comfortable "built in" headrest.
- Available in standard seat width (20") or extra wide (24")
- Available colours: Blue Ridge (BR), Jade (J) and Rosewood (R).

A41-D574-(Colour Code) . . . . . Standard Width . . . . . Each  
 A41-D574EW-(Colour Code) . . . . . Extra Wide . . . . . Each



### WINCO® 535 CAREMOR® RECLINER WITH LIQUICELL TECHNOLOGY

The 535 Caremor Recliner features infinite positions with the web suspension back and seat. The back operates using a gas spring. The release handle can be relocated underneath the left armrest. It has such standard features as convenient recessed grip areas, 5" total lock casters, non-porous blow molded tray with removable locking function and heavy gauge expanded vinyl. Colour Options: Blueridge, Royal Blue Grey, Green, Mauve, Taupe, Moss Green

509-535 . . . . . Each



### WINCO® 527 THREE POSITION RECLINER WITH DROP-ARM

The 527 Convalescent Recliner features blow-molded vinyl covered side panels with convenient recessed grip areas and pocket. It has such standard features as 5" total lock casters, non-porous blow molded tray with removable locking function, comfortable ergonomic build-in headrest and heavy gauge expanded vinyl. Colour Options: Blueridge, Burgundy, Grey, Mauve, Taupe, Moss Green

509-527 . . . . . Each



## BRODA® AUTO LOCKING GLIDER CHAIR

BRODA's Auto Locking Glider provides exceptional comfort and safety. The Auto Locking Glider provides a gentle gliding action to help create a calming, soothing effect for agitated individuals. The gliding action is excellent for Alzheimer's and Dementia as it helps reduce the need for wandering.

BRODA's unique Comfort Tension Seating® system provides pressure redistribution and air flow for increased sitting comfort and support.

Designed for Fall Prevention, the Auto Locking Glider begins to glide once the resident sits in the chair, and stops gliding once the resident begins to vacate the chair. It also features a caregiver lock to prevent unwanted gliding.

It is ideal for long-term care, assisted living and acute-care facilities for:

- Residents requiring Low to Moderate levels of nursing care.
- Alzheimer's and dementia residents unit within a facility
- Ambulatory residents who can benefit from mild exercise.
- Great for patient rooms, and lounge areas within a facility

Available in 5 different colours (Blackberry w/grey frame, Burgundy w/grey frame, Charcoal w/grey frame, Sapphire Blue w/grey frame and meteor w/ brown frame) and different seat widths and back heights.

276-100-10AL	..... 18" Seat Width & 26" Back	..... Each
276-100-15AL	..... 20" Seat Width & 26" Back	..... Each
276-100-20AL	..... 20" Seat Width & 32" Back	..... Each



## BRODA® LT TILT CHAIR

The LT helps you to provide the best quality of care for patients with removable vinyl padding that is easy to clean. Caster configuration allows for smooth tracking and ease of steering in tight corners. Small footprint compared to traditional recliners, is great in small spaces. Durable frame and components allow for overall ease of use.

### STANDARD FEATURES

- BRODA Comfort Tension Seating®
- 19" seat width, 32" back height, 20" seat height
- 28° infinitely adjustable seat tilt
- Infinitely adjustable elevating legrest
- Soft, durable polyurethane armrests

276-LT-TILT	..... 18" Width, Fixed Arms	..... Each
-------------	-----------------------------	------------



## CAREFOAM® TILT ROLLABOUT CHAIR

The gas-cylinder-assisted tilt system of this is infinitely adjustable, enabling the resident to be positioned from upright to a comfortable recline. The large wheels that are suitable for indoor and outdoor use. The Tilt Rollabout is especially recommended for those residents needing frequent position changes to reduce the possibility of tissue breakdown. The exclusive sling-back design provides lumbar support and thoracic stability, and the seating area is lined with memory foam. The adjustable leg rest is ideal for residents with contractures or for those who need maximum lower limb elevation. Available with solid blue or tapestry multi outer cover wool, blue fleece or beige fleece interior cover.



B56-CF510 ..... 18" Seat Width ..... Each  
 B56-CF520 ..... 20" Seat Width ..... Each

## CAREFOAM® ROLLABOUT CHAIR

The exclusive sling-back design accommodates the physical changes of aging while providing lumbar support and thoracic stability. This design also keeps the resident's head upright, eliminating the "looking at the ceiling" position so common with other chairs. The anti-thrust seat effectively prevents sliding out of the chair. The infinitely adjustable leg rest is ideal for residents with contractures or for those who need maximum lower limb elevation. All-foam construction and lambswool comfort cover eliminate hard edges, preventing bruising and skin tears. Available with solid blue or tapestry multi outer cover wool, blue fleece or beige fleece interior cover.



B56-CF200 ..... 18" Seat Width ..... Each  
 B56-CF220 ..... 20" Seat Width ..... Each

## CAREFOAM® TWO PIECE ISS CHAIR & OTTOMAN

The Carefoam Two-Piece Plus ISS chair features the Internal Safety System (ISS) for maximum resident protection. The ISS is a discreet, padded 3-point positioner that will ensure the safety and comfort of your resident in the most challenging situations. These chairs are appropriate tools for the care of active Huntington's clients or those with other forms of ataxia, for residents presenting with combative dementia, or for any situation where the primary concern is the safety of the resident and others. Features a detachable mobile ottoman base for ease of transport. Available with solid blue or tapestry multi outer cover wool, blue fleece or beige fleece interior cover.



B56-CF1100 ..... 18" Seat Width ..... Each  
 B56-CF1300 ..... 20" Seat Width ..... Each

## LA-Z-BOY® SLEEP CHAIR

This chair folds out and reclines in to a flat bed. Comes standard with four, tapered, wood legs or choose the optional four caster (as pictured) option for easy mobility from room to room. Chair comes standard with black urethane arm caps or you can choose the optional wood arm caps in a variety of finishes. 30" wide sleep platform.

SPFOD2800 ..... Each

\*\*Contact your Medical Mart representative for fabric options.

 **SPECIALTY**  
FURNITURE  
www.specialtyfurniture.ca



## STANDARD WHEELCHAIR WITH FIXED OR DESK ARMS

The Standard Weight Wheelchair features removable desk length (10") padded armrests or fixed full length arm rests, and swing-away footrests, which swing far out of the way for entering and exiting the chair. It also comes standard with aluminum handrims, solid tires, black upholstery, and an epoxy coated carbon steel frame.

320-MDS806150EE..... 18" Width, Fixed Arms ..... Each

320-MDS806250EE..... 18" Width, Removable Desk Arms..... Each

A19-1091..... 22" X-Wide Width, Removable Desk Arms..... Each



## TRANSPORT CHAIR

Lightweight and compact, the Transport Chair is lighter than conventional wheelchairs. Features include fold-down back, black nylon upholstery, padded armrests, 8" mag-style wheels with dual locking rear wheels, swing-away removable footrests, and standard safety belt. Easy to fold and convenient to store. 19" Wide.

320-MDS808200ARE..... Each





## PNEUMATIC HEIGHT ADJUSTABLE STOOLS

- 15" (38.1 cm) wide, 4" (10.2 cm) thick poly-foam seat cushion
- Pneumatic or spin-lift height adjustment; 17" – 21.25" (43.2 cm – 54 cm)
- Five-leg, 21" (53.3 cm) black ABS plastic composite base
- 2" (5.1 cm) black plastic dual-wheel nylon casters
- Optional adjustable contoured backrest with black plastic shroud
- 55 standard upholstery colors
- CAL-133 fire resistant upholstery available for 22500 models
- Weight capacity: 250 lb (113 kg)
- Standard 5 year warranty



## ADJUSTABLE HEIGHT STOOL, 5 CASTERS

962-22500-21	Clamshell	Each
962-22500-22	Gun Metal	Each
962-22500-24	Black	Each
962-22500-28	Cocoa	Each
962-22500-31	Deep Sapphire	Each
962-22500-32	Blue Fog	Each

Exam Stools are available in 12 standard upholstery colours. The colours above are the 6 most popular.

## ADJUSTABLE HEIGHT STOOL WITH BACKREST, 5 CASTERS

962-22500B-21	Clamshell	Each
962-22500B-22	Gun Metal	Each
962-22500B-24	Black	Each
962-22500B-28	Cocoa	Each
962-22500B-31	Deep Sapphire	Each
962-22500B-32	Blue Fog	Each



Clamshell



Gun Metal



Black



Cocoa



Deep Sapphire



Blue Fog

## RITTER® ADJUSTABLE HEIGHT STOOL, 5 CASTER

These stools offers a black, plastic composite five-leg caster base and hood, and a single lever height control. The seat height adjusts pneumatically. Black vinyl upholstery is stocked for quick delivery, other colours available with lead time.

550-273-001-(Colour Code) . . . . . Stool with Back . . . . . Each  
550-272-001-(Colour Code) . . . . . Stool without Back . . . . . Each

## RITTER® ADJUSTABLE HEIGHT STOOL, SELF LOCKING CASTERS

- Self locking casters lock when weight is applied
- Pneumatically adjusts from 18" to 24.5" in height
- Efficient air-lift design saves time and hassle
- Available with & without back
- Available in 6 standard colours (above)

550-273-005-(Colour Code) . . . . . Stool with Back . . . . . Each  
550-272-005-(Colour Code) . . . . . Stool without Back . . . . . Each

## RITTER® ADJUSTABLE HEIGHT STOOL, NON-MOVEMENT GLIDES

- Non-movement glides instead of casters
- Pneumatically adjusts from 18" to 24.5" in height
- Efficient air-lift design saves time and hassle
- Available with & without back
- Available in 6 standard colours (above)

550-273-003-(Colour Code) . . . . . Stool with Back . . . . . Each  
550-272-003-(Colour Code) . . . . . Stool without Back . . . . . Each

## REVOLVING STOOL

- Smooth screw height adjustment
- Solid steel screw shaft
- 1" chrome plated steel tubular construction
- 14" seat diameter
- 4" thick padded black vinyl seat
- Four hooded casters
- 250 lb. weight capacity

162-116-455 . . . . . Each

## STEP ON STOOLS

162-116-460 . . . . . Step On Stool . . . . . Each  
162-116-470 . . . . . Step On Stool With Handrail . . . . . Each

