

## CURAD® PAPER SURGICAL ADHESIVE TAPE

- CURAD® paper tape is a gentle, breathable tape for sensitive skin
- Recommended especially for pediatric and geriatric patients
- Flexible and conformable, hypoallergenic
- Adheres well to damp skin
- Not made with natural rubber latex

320-NON270012	..... 1/2" x 10 yard	..... 24/Box
320-NON270001	..... 1" x 10 yard	..... 12/Box
320-NON270002	..... 2" x 10 yard	..... 6/Box
320-NON270003	..... 3" x 10 yard	..... 4/Box



## CURAD® TRANSPARENT SURGICAL ADHESIVE TAPE

- CURAD® transparent tape is a breathable perforated plastic tape that permits skin examination without tape removal
- Provides a clean and easy bilateral tear
- Excellent tape to hold facial dressings or for IV sites and tubing retention
- Not made with natural rubber latex

320-NON270212	..... 1/2" x 10 yard	..... 24/Box
320-NON270201	..... 1" x 10 yard	..... 12/Box
320-NON270202	..... 2" x 10 yard	..... 6/Box



## CURAD® CLOTH SILK SURGICAL ADHESIVE TAPE

- CURAD® cloth silk-like tape is water repellent and easy to tear
- With its strong woven construction it's ideal for affixing dressings
- Not made with natural rubber latex
- Hypoallergenic

320-NON270112	..... 1/2" x 10 yard	..... 24/Box
320-NON270101	..... 1" x 10 yard	..... 12/Box
320-NON270102	..... 2" x 10 yard	..... 6/Box



### 3M MICROPORE™ PAPER TAPE

A latex-free, hypoallergenic paper tape that is gentle to the skin yet adheres well and leaves minimal adhesive residue upon removal. An economical, general purpose, breathable surgical tape

355-1530-1/2	1/2" x 10yd	24/Box
355-1530-1	1" x 10yd	12/Box
355-1530-2	2" x 10yd	6/Box



### 3M TRANSPORE™ TRANSPARENT TAPE

A latex-free, hypoallergenic, transparent and perforated plastic tape that offers strong adhesion. Transpore tape is easy to tear into very thin strips, and is easy to use even with gloves.

355-1527-1/2	1/2" x 10yd	24/Box
355-1527-1	1" x 10yd	12/Box
355-1527-2	2" x 10yd	6/Box



### 3M BLENDERM™ TRANSPARENT TAPE

A latex-free, hypoallergenic, occlusive and transparent plastic tape. This flexible medical tape is designed to provide a barrier from external fluids and contamination.

355-1525-1/2	1/2" x 5yd	24/Box
355-1525-1	1" x 5yd	12/Box
355-1525-2	2" x 5yd	6/Box



### 3M KIND REMOVAL SILICONE TAPE

Softer and more conformable than traditional adhesives, works differently with skin and hair. Can be removed with minimal stripping of skin cells or pulling of hair. Constant adhesive strength for as long as it's left in place, unlike traditional tape that increases adhesion over time.

355-2770-1	1" x 5.5 Yard	12/Box
355-2770-2	2" x 5.5 Yard	6/Box



### ADHESIVE CLOTH TAPE

366-202001	1" x 10yd	Box/12
366-202002	2" x 10yd	Box/6



### ADHESIVE REMOVER

Removes tape and sticky residue

#### Smith & Nephew Remove

368-403379	Remove® 220mL	Each
368-403120	Remove® Prep Pad	50/Box

#### Generic

320-MDS090855	Prep Pad	100/Box
---------------	----------	---------



## MEDFIX® DRESSING RETENTION TAPE

- For flexibility and customized sizing.
- Low sensitivity adhesive, gentle for the patient
- Medfix™ has a printed s-c release liner
- Water resistant
- Indications: to secure primary dressings and to secure gastrostomy tubes and other feeding tubes
- Contraindicated as a primary dressing
- May be left in place for up to 7 days
- Dressing change frequency will depend on amount of drainage

320-MSC4002	.....	5cm x 10 Meter	.....	Each
320-MSC4004	.....	10cm x 10 Meter	.....	Each
320-MSC4006	.....	15cm x 10 Meter	.....	Each



## MEDFIX® EZ DRESSING RETENTION TAPE

- Low sensitivity adhesive, gentle for the patient
- Medfix™ EZ is linerless with 2" perforations
- Water resistant
- Indications: to secure primary dressings and to secure gastrostomy tubes and other feeding tubes
- Contraindicated as a primary dressing
- May be left in place for up to 7 days
- Dressing change frequency will depend on amount of drainage

320-MSC4102	.....	5cm x 10 Meter	.....	Each
320-MSC4104	.....	10cm x 10 Meter	.....	Each
320-MSC4106	.....	15cm x 10 Meter	.....	Each



## GENTAC® SILICONE DRESSING FIXATION TAPE

- Soft silicone adhesive ideal for fixation on fragile skin
- Gentle for patient
- Cut to size needed; easy to apply
- Waterproof
- Indications: to secure primary or secondary dressings and to secure gastrostomy tubes and other feeding tubes
- Contraindicated as a primary dressing
- May be left in place for up to 7 days; dressing change frequency will depend on amount of drainage

320-MSC1583	.....	2cm x 3 Meter	.....	12/Box
320-MSC1585	.....	5cm x 3 Meter	.....	6/Box



## MEFIX® TAPE

Mefix tape by Molnlycke can be used in a wide variety of situations where fixation is required, including fixation of dressings, swabs, catheters, and tubes. The porous structure of the non-woven fabric ensures air and water vapor permeability. The fabric is flexible, conforms well to body contours and allows body movements when in place. The fabric is non-absorbent and provides a barrier for blood and fluid strike-through.

173-310250	2.5cm x 10m	Each
173-310500	5cm x 10m	Each
173-311000	10cm x 10m	Each
173-311500	15cm x 10m	Each
173-312000	20cm x 10m	Each



## HYPAFIX® TAPE

A cross-elastic non-woven tape that provides light compression. Excellent for areas requiring flexibility and freedom from constraint while providing superior holding power. Air and exudate permeable. Hypoallergenic.

127-7144300	2.5cm x 10m	Each
127-7144301	5cm x 10m	Each
127-7144302	10cm x 10m	Each
127-7144303	15cm x 10m	Each
127-7144304	20cm x 10m	Each



## 3M SOFT CLOTH ON A LINER

The same soft, latex-free, hypoallergenic, conformable tape as standard Medipore tape. Gentle enough for at-risk skin yet strong enough for critical securement, this tape features an s-slit liner with a printed grid and can be easily cut to a desired size.

355-2762	5cm x 9.1m	Each
355-2764	10cm x 9.1m	Each



## 3M MEDIPORE-H™ SOFT CLOTH TAPE

Soft, gentle and conformable tape, offering enhanced gentleness and stronger adhesion for challenging tape applications while being easy to use and easy on the skin. Ideal for anchoring post-surgical dressings where swelling is expected.

355-2861	2.5cm x 9.1m	2/Package
355-2862	5cm x 9.1m	Each
355-2863	7.5cm x 9.1m	Each
355-2864	10cm x 9.1m	Each
355-2866	15cm x 9.1m	Each
355-2868	20cm x 9.1m	Each



## MEPITAC® ADHESIVE TAPE

Mepitac is highly conformable and easy to apply and remove. The Safetac layer assures secure fixation and gentle removal from skin and bodily hair. If inspection or repositioning is necessary, Mepitac can be lifted and re-applied. Mepitac is an ideal fixation solution for patients with fragile and sensitive skin. Mepitac is suitable for patients requiring repeated application and removal of tape over the same area.

173-298300	2cm x 3m	Each
173-298400	4cm x 1.5m	Each



**PLASTIC STRIP ADHESIVE BANDAGES  
3/4" X 3"**



211-018-600 ..... 100/Box

**PLASTIC STRIP ADHESIVE BANDAGES  
1" X 3"**



320-NON25600 ..... 100/Box

**PLASTIC SPOT BANDAGE  
7/8" ROUND**



211-P150150..... 100/Box

**FABRIC ADHESIVE STRIP BANDAGES**



3/4" x 3"

320-NON25650 ..... 100/Box

1" x 3"

320-NON25660 ..... 100/Box

**FABRIC FINGERTIP BANDAGE**



1.75" x 2"

320-NON25513 ..... 100/Box

**FABRIC KNUCKLE BANDAGE**



1.5" x 3"

320-NON25510..... 100/Box

**TRIANGULAR BANDAGE**



211-PM3-600..... Each

**STERION® STRIP SKIN CLOSURES**



294-75-1-040PCS...1/8" x 3" .....250/Box

294-75-1-041PCS...1/4" x 3" .....150/Box

294-75-1-042PCS...1/4" x 1-1/2" ...300/Box

294-75-1-046PCS...1/4" x 4" .....150/Box

294-75-1-047PCS...1/2" x 4" .....150/Box

**COVERPLAST® CLASSIC FABRIC DRESSING STRIPS**



127-7115104 .....4cm x 5m.....Each

127-7115103 .....6cm x 5m.....Each

127-7115102 .....8cm x 5m.....Each

## BORDERED GAUZE DRESSING

- Deluxe soaker non-adherent pad
- Non-woven adhesive border with water resistant backing
- Ideal for incision sites
- Indications: Pressure ulcers, Partial and full-thickness wounds, Incision sites

320-MSC3222	5cm x 5cm with 2.5cm x 2.5cm Pad	15/Box
320-MSC3244	10cm x 10.5cm with 6cm x 6cm Pad	15/Box
320-MSC3245	10cm x 12.5cm with 5cm x 6cm Pad	15/Box
320-MSC3248	10cm x 20cm with 5cm x 15cm Pad	15/Box
320-MSC32410	10cm x 25cm with 5cm x 20cm Pad	15/Box
320-MSC3266	15cm x 15cm with 10cm x 10cm Pad	15/Box



## SURESITE 123+ PAD® TRANSPARENT FILM DRESSINGS

- Easy delivery of an all-in-one cover dressing with absorption
- Traditional moisture vapor transmission rate (MVTR)
- Conformable; won't stick to itself
- Microporous technology
- Permits continuous observation
- Indications: partial-thickness wounds, full-thickness wounds (secondary dressing), peripheral and central IV lines, skin tears, lacerations and abrasions and to help prevent skin breakdown caused by friction to epidermis
- Contraindications: contraindicated as a primary dressing on wounds with moderate to heavy drainage

320-MSC2603	6cm x 7cm	100/Box
320-MSC2605	10cm x 12cm	50/Box
320-MSC2610	8.9cm x 25.4cm	25/Box



## SUREFILM® TRANSPARENT FILM ROLLS

- SureSite® transparent roll allows for customization of size and shape of transparent film
- Traditional moisture vapor transmission rate (MVTR)
- Conformable; won't stick to itself
- Microporous technology
- Permits continuous observation
- Indications: partial-thickness wounds, full-thickness wounds (secondary dressing), peripheral and central IV lines, skin tears, lacerations and abrasions and to help prevent skin breakdown caused by friction to epidermis
- Contraindications: contraindicated as a primary dressing on wounds with moderate to heavy drainage

320-MSC2402	5cm x 10 meter	Each
320-MSC2404	10cm x 10 meter	Each



### NON-ADHERENT PADS (TELFA TYPE)

320-NON25700	2"x3"	100/Box
510-6017	2"x3" With adhesive	100/Box
320-NON25710	3"x4"	100/Box
320-NON25720	3"x8"	50/Box



### TELFA PLUS® ADHESIVE ISLAND NON-ADHERENT PADS

510-7540	4" x 5"	25/Box
510-7541	4" x 8"	25/Box



### 3M TEGADERM™ + PAD DRESSINGS

Waterproof, sterile, all-in-one dressing with non-adherent absorbent pad. Ideal for surgical wounds and small incisions or wounds.

355-3582	5cm x 7cm	50/Box
355-3584	6cm x 10cm	50/Box
355-3586	9cm x 10cm	25/Box
355-3587	8cm x 11cm Oval	25/Box
355-3588	15cm x 15cm	25/Box
355-3589	9cm x 15cm	25/Box
355-3590	9cm x 20cm	25/Box
355-3591	9cm x 25cm	25/Box



### 3M MEDIPORE™ + PAD DRESSINGS

Soft, conformable, all-in-one sterile dressing with non-adherent absorbent pad. Ideal for surgical wounds and small incisions or wounds.

355-3562	5cm x 7cm	50/Box
355-3564	6cm x 10cm	50/Box
355-3566	8cm x 10cm	25/Box
355-3569	9cm x 15cm	25/Box
355-3570	9cm x 20cm	25/Box
355-3571	9cm x 25cm	25/Box
355-3573	9cm x 35cm	25/Box



## PACKING STRIPS, STERILE, PLAIN & WITH IODOFORM

### PLAIN, 1 per Bottle

819-3401 .....	1/4" x 5yds
819-3402 .....	1/2" x 5yds
819-3403 .....	1" x 5yds
819-3404 .....	2" x 5yds

### WITH IODOFORM, 1 per Bottle

819-3411 .....	1/4" x 5yds
819-3412 .....	1/2" x 5yds
819-3413 .....	1" x 5yds
819-3414 .....	2" x 5yds



## COVIDIEN CURITY AMD PACKING STRIP

Curity™ AMD packing strips are impregnated with PHMB (Polyhexamethylene Biguanide), a powerful yet safe antiseptic that has a broad range effectiveness against gram positive and gram negative micro-organisms including some multi-drug resistant strains such as MRSA. This exciting new product provides healthcare facilities a low cost, prophylactic method of treatment without changing nursing protocol.

510-7831AMD .....	1/4" x 36"	.....	10/Box
510-7832AMD .....	1/2" x 36"	.....	10/Box
510-7833AMD .....	1" x 36"	.....	10/Box



## COVIDIEN KERLIX AMD ANTI-MICROBIAL BANDAGE, STERILE

Kerlix™ AMD rolls are impregnated with PHMB (Polyhexamethylene Biguanide), a powerful yet safe antiseptic that has a broad range effectiveness against gram positive and gram negative micro-organisms including some multi-drug resistant strains such as MRSA. This exciting new product provides healthcare facilities with a low cost, prophylactic method of treatment without changing nursing protocol. In addition to its antimicrobial capabilities, Kerlix™ AMD rolls provide fast-wicking action, superior aeration, and excellent absorbency.

510-3332 .....	4.5" x 4.1 yards	.....	Each
----------------	------------------	-------	------



## COVIDIEN KERLIX AMD ANTI-MICROBIAL SUPER GAUZE/SPONGE, STERILE

Kerlix™ AMD antimicrobial dressings super sponges are impregnated with PHMB (Polyhexamethylene Biguanide), a powerful yet safe antiseptic that has a broad range effectiveness against gram positive and gram negative micro-organisms including some multi-drug resistant strains such as MRSA. This exciting new product provides healthcare facilities a low cost, prophylactic method of treatment without changing nursing protocol. In addition to its antimicrobial capabilities the Kerlix™ AMD antimicrobial dressings super sponges are a superior primary and secondary dressing designed to absorb and transfer exudate away from the wound.

510-6665 .....	15cm x 17cm	.....	5/Pkg
----------------	-------------	-------	-------





## WOVEN GAUZE SPONGES, NON-STERILE

### 8-Ply Non-Sterile

211-052113	2" X 2"	Pkg/200
211-052114	3" X 3"	Pkg/200
211-052110	4" X 4"	Pkg/200

### 12-Ply Non-Sterile

211-052113	2" X 2"	Pkg/200
211-052122	3" X 3"	Pkg/200
211-052123	4" X 4"	Pkg/200



## WOVEN GAUZE SPONGES, STERILE 2'S

320-NON21420	2" X 2"	50 Pkg/2/Box
320-NON21422	3" X 3"	40 Pkg/2/Box
320-NON21448	4" X 4"	50 Pkg/2/Box



## NON-WOVEN GAUZE SPONGES 4-PLY, NON-STERILE

211-94122	2" X 2"	Pkg/200
211-94133	3" X 3"	Pkg/200
211-94144	4" X 4"	Pkg/200



## NON-WOVEN GAUZE SPONGES, 4 PLY, STERILE 2'S

211-A2240	2" X 2"	Box/50 Pkg of 2
211-A3340	3" X 3"	Box/50 Pkg of 2
211-A4440	4" X 4"	Box/50 Pkg of 2
211-A8840	4" X 8"	Box/50 Pkg of 2



## DRAIN/TRACH SPONGE, NON-WOVEN, STERILE 2'S

211-A9210	4" x 4", 6 ply	25 x 2/Box
-----------	----------------	------------



**CONFORM STRETCH GAUZE, NON-STERILE**

211-A205	2" x 5 yds	12/Box
211-A305	3" x 5 yds	12/Box
211-A405	4" x 5 yds	12/Box
211-A605	6" x 5 yds	6/Box



**CONFORM STRETCH GAUZE, STERILE**

211-S134CA/50	2" x 5 yds	12/Box
211-S134CA/75	3" x 5 yds	12/Box
211-S134CA/100	4" x 5 yds	12/Box
211-S134CA/150	6" x 5 yds	6/Box



**ELASTIC TENSOR BANDAGES**

211-A620	2" x 5yds	Each or 10/Box
211-A621	3" x 5yds	Each or 10/Box
211-A622	4" x 5yds	Each or 10/Box
211-A623	6" x 5yds	Each or 10/Box



**ABDOMINAL PADS**

211-A7059	Sterile, 5" x 9"	20/Tray
320-NON21451	Non-Sterile, 5" x 9"	576/Case
211-A7081	Sterile, 8" x 10"	24/Tray



**COMBINE ROLL**

238-A7182	8"x20yd	Each
-----------	---------	------



**ELASTIC NET BANDAGE**

320-NONNET01	Size 1	25yds/Box
320-NONNET02	Size 2	25yds/Box
320-NONNET03	Size 3	25yds/Box
320-NONNET04	Size 4	25yds/Box
320-NONNET05	Size 5	25yds/Box
320-NONNET06	Size 6	25yds/Box
320-NONNET07	Size 7	25yds/Box
320-NONNET08	Size 8	25yds/Box
320-NONNET09	Size 9	25yds/Box
320-NONNET10	Size 10	25yds/Box



## MEDIGRIP® ELASTIC TUBULAR BANDAGE

- Provides excellent support for joints
- Easy to apply and reapply
- Wide range of applications
- Good for securing dressings
- Can be used as mild compression when doubled
- Indications: edema, treatment of chronic venous insufficiency, dislocations, sprains and hypertrophic scarring
- Contraindicated for patients with a known sensitivity to latex
- Dressing change frequency will depend on amount of drainage
- Each roll is 11 yards (10 meters) in length
- Size A measures 1-3/4"W (4.5 cm), for infant feet and arms; Size B measures 2-1/2"W (6.3 cm), for small hands and limbs; Size C measures 2-5/8"W (6.8 cm), for adult hands, arms or legs; Size D measures 3"W (7.5 cm), for large arms or legs; Size E measures 3-1/2"W (8.75 cm), for legs or small thighs; Size F measures 4"W (10 cm), for large knees or thighs; Size G measures 4-3/4"W (12 cm), for large thighs.



320-MSC9500.....	Size A.....	4.5cm.....	10m Roll/Box
320-MSC9501.....	Size B.....	6.3cm.....	10m Roll/Box
320-MSC9502.....	Size C.....	6.8cm.....	10m Roll/Box
320-MSC9503.....	Size D.....	7.5cm.....	10m Roll/Box
320-MSC9504.....	Size E.....	8.75cm.....	10m Roll/Box
320-MSC9505.....	Size F.....	10cm.....	10m Roll/Box
320-MSC9506.....	Size G.....	12cm.....	10m Roll/Box

## COFLEX® TLC COMPRESSION SYSTEM

- Advanced two layer compression system provides support and compression to feet and legs
- Layer 1 is a soft foam with an absorbent coating (blue contact layer) that wicks away moisture and helps control odor
- Layer 2 is a short stretch compression layer that provides adaptive compression levels and can be easily torn by hand
- Both layers are not made with natural rubber latex and bond together when applied which keeps the system in place
- A nylon stocking is included to apply over the completed dressing for patient comfort and ease of movement under clothes and on bed sheets
- Comfortable and lightweight for increased patient compliance



320-AND7800.....	Coflex TLC Kit.....	Each
320-AND7800.....	Coflex TLC Kit.....	8 Kits/Case
320-AND7800XL.....	Coflex TLC XL Kit.....	Each
320-AND7800XL.....	Coflex TLC XL Kit.....	8 Kits/Case
320-AND7802.....	Coflex TLC Lite Kit.....	Each
320-AND7802.....	Coflex TLC Lite Kit.....	8 Kits/Case

### THREEFLEX® BANDAGE SYSTEM

- For lighter compression or for mixed etiology
- Effective appropriate compression
- Extended wear time
- Absorbs drainage
- Efficient packaging
- Contraindicated for patients with severe arterial disease
- Multi-layer compression bandages may be left in place for up to 7 days
- Dressing change frequency will depend on amount of drainage
- Kit includes: Padding (4" x 3.8 yd, 10.2 cm x 3.5 m) unstretched; Light Conforming (4" x 3 yd, 10.2 cm x 2.7 m) unstretched; Cohesive (4" x 3.2 yd, 10.2 cm x 2.9 m) unstretched; Medi-Strips



320-MSC4300	.....	Threeflex System	.....	Each
320-MSC4300	.....	Threeflex System	.....	8 Kits/Case

### FOURFLEX® BANDAGE SYSTEM

- For the treatment of Chronic Venous Insufficiency
- Effective appropriate compression
- Extended wear time
- Absorbs drainage
- Efficient packaging
- Contraindicated for patients with severe arterial disease
- Multi-layer compression bandages may be left in place for up to 7 days
- Dressing change frequency will depend on amount of drainage
- Kit includes: Padding (4" x 3.8 yd, 10.2 cm x 3.5 m) unstretched; Light Conforming (4" x 3 yd, 10.2 cm x 2.7 m) unstretched; Compression (4" x 5 yd, 10.2 cm x 4.6 m) unstretched; Cohesive (4" x 3.2 yd, 10.2 cm x 2.9 m) unstretched; Medi-Strips



320-MSC4400	.....	Fouflex System	.....	Each
320-MSC4400	.....	Fourflex System	.....	8 Kits/Case
320-MSC4400XL	.....	Fouflex System	.....	Each
320-MSC4400XL	.....	Fourflex System	.....	8 Kits/Case

### STOCKINETTE TUBULAR 100% COTTON BANDAGE

211-459220	2"	25yds/Box
211-459240	3"	25yds/Box
211-459250	4"	25yds/Box
211-459270	6"	25yds/Box



### COMPRILAN® COMPRESSION BANDAGE

127-7718700	4cm x 5m	Each
127-0102600	6cm x 5m	Each
127-0102700	8cm x 5m	Each
127-0102800	10cm x 5m	Each
127-0102900	12cm x 5m	Each



### EASIFIX® COHESIVE BANDAGE

127-7143610	6cm x 4m	10/Box
127-7143611	8cm x 4m	10/Box
127-7143612	10cm x 4m	10/Box
127-7143613	12cm x 4m	10/Box
127-7143614	4cm x 20m	Each
127-7143615	6cm x 20m	Each
127-7143616	8cm x 20m	Each
127-7143617	10cm x 20m	Each
127-7143618	12cm x 20m	Each



### TUBIGRIP® TUBULAR BANDAGE

173-1436	Size B	10m/Box
173-1443	Size C	10m/Box
173-1437	Size D	10m/Box
173-1434	Size E	10m/Box
173-1438	Size F	10m/Box
173-1439	Size G	10m/Box
173-1440	Size J	10m/Box



### TUBIFAST® TUBULAR BANDAGE

173-2434	3.5cm x 10m, Small Limbs	Each
173-2436	5cm x 10m, Small & Medium Limbs	Each
173-2438	7.5cm x 10m, Large Limbs	Each
173-2440	10.75cm x 10m, X-Large Limbs	Each



**DRESSING TRAY, STERILE**

- |                              |                          |
|------------------------------|--------------------------|
| 3 Compartment Tray Contains: | 6 ea Cotton Balls        |
| 3 ea Plastic Forceps         | 1 ea Waste Bag with Ties |
| 2 ea CSR Towels              | 6 ea 3" x 3" Gauze       |
- 211-D101108.....Each  
 211-D101108.....100/Case



**DRESSING TRAY, STERILE**

- |                                  |                                   |
|----------------------------------|-----------------------------------|
| 3 Compartment Tray Contains:     | 2 Towels                          |
| 2 Metal thumb forceps            | 1 Waste bag with tie              |
| 6 Gauze sponges 5cm x 5cm, 8 ply | 2 Gauze sponges 10cm x 10cm 8 ply |
| 1 Plastic transfer forceps       |                                   |
- 320-M101161A.....20/Case



**SUTURE REMOVAL TRAY, STERILE**

- 320-MDS707555 .....50/Case



**TRACHEOSTOMY TRAY, STERILE**

- Contains:
- |                            |                               |
|----------------------------|-------------------------------|
| - 1 ea 500mL Tray          | - 2 ea Cotton Tip Applicators |
| - 1 ea Cleaning Brush      | - 1 ea Twill Tape             |
| - 2 ea 4 x 4 Gauze Sponges | - 1 ea Tracheostomy Dressing, |
| - 2 ea Pipe Cleaners       | - 1 ea Waterproof Drape       |
| - 1 Pair of Vinyl Gloves   | - 1 ea Pop-up Basin           |
- 189-AS860.....Each



**LACERATION TRAY, STERILE**

- Contains:
- |   |                                  |
|---|----------------------------------|
| - 2-Compartment Tray,                     | - 1 ea Curved Mosquito Forceps   |
| - 5 ea 8-ply 2" x 2" Gauze                | - 1 ea Webster Needle Holder     |
| - 10 ea 8-ply 4" x 4" Gauze               | - 2 ea 2 oz (59mL) Medicine Cups |
| - 2 ea Blue Blotting Towels (19.5" x 21") | - 1 ea 10mL Syringe              |
| - 1 ea Fenestrated Drape                  | - 1 ea 18g x 1.5 in Needle       |
| - 1 ea Adson Thumb Forceps 1 x 2          | - 1 ea 25g x 5/8 in Needle       |
| - 1 ea Straight Iris Scissors             | - 1 ea 27g x 1/2 in Needle.      |
- 320-DYNJ03000 ..... Each



**SKIN STAPLER SET, 35 SHOT**



344-PMW35 . . . . Each

**SKIN STAPLE REMOVER KIT, STERILE**



402-716 . . . . Each

**FIRST AID KITS**

822-50424 . . . . . 6-15 Employee . . . . . Each  
 485-F852M021 . . . . . 16-200 Employee . . . . . Each



**EYE WASH BOTTLE & SOLUTION**

485-F4560701 . . . . . Eyewash Station with 1 L bottle . . . . . Each  
 485-F4501169 . . . . . Eyewash Solution, 1 Litre . . . . . Each  
 485-F4570301 . . . . . Eyewash Bottle Only, 1 L . . . . . Each



**EYE PADS, STERILE**

320-NON21600 . . . . . 1-5/8" x 2-5/8" . . . . . 50/Box



**POLYSPORIN® OINTMENT**

189-23789 . . . . . 28gr Tube . . . . . Each  
 189-23798 . . . . . 15gr Tube . . . . . Each



**CARE CADDY**

Small Rubbermaid carry caddy with handle is ideal for nursing staff to have all essential supplies for skin & wound care in to carry to the bedside for easy access

700-3154-88 . . . . . 16"w x 11" d x 6.75"h . . . . . Each

