

BD VERITOR® FLU TEST KIT

Testing for Influenza using the BD Veritor™ System provides you with a timely and accurate diagnosis, ensuring residents receive the right treatment earlier (reduction in use of antivirals), while also allowing you to quickly quarantine sick residents in your homes thus preventing the potential spread of the disease.

- Simple - Less than 50 seconds of hands-on time with an objective, digitally displayed result
- Accurate - The first Digital Immunoassay (DIA), a new category of diagnostic tests that combines advances in detection particles, optical image recognition, and interpretation algorithms to improve accuracy
- Fast - From initial sample processing to a digital result in minutes

| | | |
|------------|---------------------------|--------|
| 308-256045 | Veritor Flu A+B Test Kits | Box/30 |
| 308-256055 | Veritor System Reader | Each |



URISYS 1100® URINE ANALYZER

With simple, 1-button operation, the Urisys 1100 Urine Analyzer makes urine testing fast, easy and accurate. Urinalysis can be a valuable diagnostic tool and a time-intensive one. With the push of a button, it evaluates the test strip and prints a patient report.

- Automatic result printout ensures consistent, accurate results
- Simple cleaning minimizes maintenance needs
- Flexible strip options: Chemstrip® 10 MD*, Chemstrip 7 or Chemstrip 5 OB test strips
- Error messages alert user to issues that could affect accuracy of results

| | | |
|-----------------|--|--------|
| 711-03617548001 | Urine strip analyzer with chemstrip technology | Each |
| 711-06431321001 | Urisys 1100 printer paper roll | Box/20 |
| 711-11248693001 | Urisys 1100 keyboard | Each |
| 711-11379194263 | Urisys 1100 control test | Box/50 |



CHEMSTRIP™ URINE TEST STRIPS

| | | |
|-----------------|--------------|------------|
| 711 1203479 | Chemstrip 10 | Bottle/100 |
| 711 10398411119 | Chemstrip 9 | Bottle/100 |
| 711 11008552119 | Chemstrip 7 | Bottle/100 |
| 711 11893467119 | Chemstrip 5 | Bottle/100 |
| 711 10200743119 | Chemstrip GP | Bottle/100 |



SIEMENS URINE TEST STRIPS

| | | |
|-----------|-----------------|------------|
| 553-2300A | Multistix 10 SG | Bottle/100 |
| 553-2304A | Multistix 8 SG | Bottle/100 |
| 553-2309 | Multistix 5 | Bottle/100 |
| 553-2810A | Labstix | Bottle/100 |
| 553-2855U | Uristix | Bottle/100 |



COAGUCHEK® INR MONITOR

The smart way to test INR

The CoaguChek XS system is a dedicated solution for patients on vitamin K antagonist therapy to monitor their PT/INR values. Using a small drop (8 µl) of capillary blood, it is easy to use and it provides a precise PT/INR result usually in less than 1 minute.

- Saves caregiver time retrieving calls from lab
- No early morning 'wake-up' blood draws for resident
- Finger lance vs painful venous draw
- Improved time in therapeutic range for resident which leads to improved outcomes
- Physicians can adjust their own patient's dosages the day they are on site

| | | |
|-----------------|---|---------------|
| 711-04870026001 | CoaguCheck XS Monitor system | Each |
| 711-04625374190 | CoaguCheck PT test strip for XS Monitor | Box/6 Strips |
| 711-04625358019 | CoaguCheck PT test strip for XS Monitor | Box/24 Strips |
| 711-04625315019 | CoaguCheck PT test strip for XS Monitor | Box/48 Strips |
| 711-03506509001 | CoaguChek Lancets, yellow | Box/50 |
| 711-04696522190 | CoaguChek XS PT Control Solution 0.8ml | Box/4 |



ONETOUCH® ULTRA BLOOD GLUCOSE TESTING

Inside every OneTouch Ultra Blue Test Strip, DoubleSure Technology measures each blood sample not once, but twice. So you're not just sure of your result, you're DoubleSure.

| | | |
|-------------|-------------------------------------|---------|
| 044-021-130 | One Touch Ultra blood glucose meter | Each |
| 044-022-901 | One Touch Ultra test strips | 100/Box |
| 044-020-348 | One Touch Ultra control solution | 2/Box |



BAYER® CONTOUR® NEXT BLOOD GLUCOSE TESTING

Innovative technology; giving you results that are even closer to those obtained in a laboratory. Bayer's test strips code the meter automatically for every test, and need a small blood sample (0.6 µL). The meter comes with a large, easy-to-read display and has a 5-second result countdown.

| | | |
|----------|--------------------------|---------|
| 032-7268 | Contour Next EZ meter | Each |
| 032-7322 | Contour Next test strips | 100/Box |
| 032-7325 | Control Solution | Each |



ACCUTREND® GC CHOLESTEROL & GLUCOSE TESTING

The Accutrend GC meter is both fast and easy to use and is most suited for accurate screening, diagnosis and monitoring of both cholesterol and glucose levels.

| | | |
|-----------------|--|--------|
| 711-11418262171 | Accutrend cholesterol test strips | 25/Box |
| 711-11447475187 | Accutrend glucose test strips | 25/Box |
| 711-11284878190 | Accutrend glucose control solution | 2/Box |
| 711-11418289190 | Accutrend cholesterol control solution | 1/Box |

